



The Importance of Quality Assurance Roles in International Trade

Testing for antibiotics in food is essential for international trade. Different countries may have different regulations regarding antibiotic residues in food. By conducting testing, food producers can demonstrate compliance with these regulations, facilitating the export and import of food products across borders.

See What Our Customer Says About Us

AB Labo France

Analyser: Evidence Investigator

Kits: Antimicrobial Array I Ultra / Antimicrobial Array II Plus / Antimicrobial Array III /
10mM 4-Nitrobenzaldehyde (in DMSO) 10ml / Honey

“We have been using the Evidence Investigator daily for 10 years now. Randox Food Diagnostic’s professionalism and customer support are major assets in the bond of trust between our two companies.”



The Randox Evidence Investigator is a multi-analyte semi-quantitative drug residue screening analyser offering laboratory results comparable to that of LC-MS/MS and other confirmatory methods.

The Evidence Investigator runs Biochip Array Technology to help prevent antibiotic resistance by offering multiple screening solutions across the following industries: honey, meat, seafood, feed and milk.

Biochip Array Technology is a platform that can provide the screening of up to 48 food or feed samples, providing results for veterinary drug residues and toxins in under 3 hours.

